

Product Datasheet

MAGED1 antibody (orb783067)

Description	MAGED1 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human MAGED1 (50-200/778 aa)
Target	MAGED1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	85 kDa
Uniprot ID	Q9Y5V3
Dilution Range	WB=1:500-1000
Expiration Date	12 months from date of receipt.