



Product Datasheet

Human Kappa Light Chain antibody (orb783134)

Description	Human Kappa Light Chain antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Full length native human IgG (1-107/107 aa)
Target	Human Ig Kappa light chain
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	11 kDa
Uniprot ID	P01834
Dilution Range	WB=1:500-2000, IHC-P=1:400-800
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human
breast
carcinoma...



Western
blot: Lane
1: Human
Serum; Lane
2: H...