

Product Datasheet

ERK1 antibody (orb783151)

Description	ERK1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from mouse p44/42 MAPK around the phosphorylation site of Thr202/Tyr204 (FL(p-T)E(p-Y)VA)
Target	phospho-ERK1/2 (Thr202 + Tyr204)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	41 kDa
Uniprot ID	Q63844
Expiration Date	12 months from date of receipt.