

Product Datasheet

RNF4 antibody (orb783210)

Description	RNF4 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human RNF4.
Target	RNF4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	21 kDa
Uniprot ID	P61020
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.