

Product Datasheet

GOLGA5 antibody (orb783246)

Description	GOLGA5 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human GOLGA5 between 1-300 amino acid.
Target	GOLGA5
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	71 kDa
Uniprot ID	Q8TBA6
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.