

Product Datasheet

SULT2B1 antibody (orb783339)

Description	SULT2B1 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human SULT2B1.
Target	SULT2B1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	41 kDa
Uniprot ID	O00204
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.