

Product Datasheet

ME2 antibody (orb783348)

Description	ME2 antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ME2 (501-584/584 aa)
Target	ME2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	63 kDa
Uniprot ID	P23368
Dilution Range	WB=1:500-4000, IHC-P=1:25
Expiration Date	12 months from date of receipt.