

Product Datasheet

LC3 antibody (orb783388)

Description	LC3 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human MAP1LC3A/B
Target	LC3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	13 kDa
Uniprot ID	Q9GZQ8
Dilution Range	WB=1:500-1000
Expiration Date	12 months from date of receipt.