



Product Datasheet

CD73 antibody (orb783396)

Description	CD73 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Recombinant human CD73 protein
Target	CD73
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	63 kDa
Uniprot ID	Q61503
Dilution Range	WB=1:300-1000, IHC-P=1:200-400, Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.



Flow cytometry of A375.Primary Antibody ...



IHC-P analysis of human colon carcinoma ...