

Product Datasheet

VCP antibody (orb783397)

Description	VCP antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	Recombinant human VCP.
Target	VCP
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	97 kDa
Uniprot ID	P55072
Dilution Range	WB=1:500-2000, IHC-P=1:25, ICC=1:25, Flow-Cyt=1:25, IHC-F=1:25, IF=1:25
Expiration Date	12 months from date of receipt.