

Product Datasheet

CDK5 antibody (orb783398)

Description	CDK5 antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	Recombinant human CDK5.
Target	CDK5
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	32 kDa
Uniprot ID	Q00535
Dilution Range	WB=1:500-2000, IHC-P=1:25, ICC=1:25, Flow-Cyt=1:25, IHC-F=1:25, IF=1:25
Expiration Date	12 months from date of receipt.