


Product Datasheet

BCL10 antibody (orb783402)

Description	BCL10 antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human BCL10 (1-130/233 aa)
Target	BCL10
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	26 kDa
Uniprot ID	O95999
Dilution Range	WB=1:500-2000, IHC-P=1:400-800, Flow-Cyt=1:100, IHC-F=1:400-800
Expiration Date	12 months from date of receipt.



Flow
cytometry of
Hela.Primary
Antibody ...