

Product Datasheet

RPS6KB2 antibody (orb783404)

Description	RPS6KB2 antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human RPS6KB2 (His-tagged).
Target	RPS6KB2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1,K
Clonality	Monoclonal
MW	53 kDa
Uniprot ID	Q9UBS0
Dilution Range	WB=1:500-2000, IHC-P=1:200-500
Expiration Date	12 months from date of receipt.