

Product Datasheet

AKT2 antibody (orb783405)

Description	AKT2 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human AKT2.
Target	AKT2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	56 kDa
Uniprot ID	P31751
Dilution Range	WB=1:500-1000, IHC-P=1:50-100
Expiration Date	12 months from date of receipt.