


Product Datasheet

MMP2 antibody (orb783406)

Description	MMP2 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	Recombinant human MMP2.
Target	MMP2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG2b,k
Clonality	Monoclonal
MW	72 kDa
Uniprot ID	P08253
Dilution Range	WB=1:500-1000, IHC-P=1:50, ICC=1:50
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human
gastric
carcinom...