


Product Datasheet


HMGB1 antibody (orb783417)

Description	HMGB1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human HMGB1 (100-200/215 aa)
Target	HMGB1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	25 kDa
Uniprot ID	P09429
Dilution Range	WB=1:500-2000, IHC-P=1:50-100, IF=1:50-100
Expiration Date	12 months from date of receipt.




IHC-P
analysis of
rat brain
using HMGB1

...



IHC-P
analysis of
mouse liver
using
HMGB...



Western
blot: Lane
1: Lung
(Mouse); Lane
2: ...