



Product Datasheet

MARK2 (phospho-Ser595) antibody (orb783419)

Description	MARK2 (phospho-Ser595) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human MARK2 around the phosphorylation site of Ser595 (SR(p-S)TF)
Target	phospho-MARK2 (Ser595)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	82-90 kDa
Uniprot ID	Q7KZI7
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
Rat liver
using p-
MARK...



IHC-P
analysis of
Human
kidney
using pho...