


Product Datasheet

**GSK3 Alpha + Beta (phospho-Tyr279+Tyr216)
antibody [Out of stock] [Out of stock] (orb783420)**

Description

GSK3 Alpha + Beta (phospho-



Western blot:brain
(Mouse);Primary:Anti-

...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, WB

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human GSK-3 Beta around the phosphorylation site of Tyr216 (VS(p-Y)IC)

Target

phospho-GSK3 Alpha + Beta
(Tyr279+Tyr216)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin,
0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20°C for one
year. Avoid repeated freeze/thaw
cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

47/51 kDa

Uniprot ID

P49841

Dilution Range

WB=1:500-2000, Flow-Cyt=1µg /test,
IHC-P=1:100-500, IHC-F=1:100-500,
IF=1:100-500, ELISA=1:5000-10000

Expiration Date

12 months from date of receipt.