



Product Datasheet

MBP antibody (orb783422)


Description	MBP antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human MBP (241-304/304 aa)
Target	MBP
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	33 kDa
Uniprot ID	P02686
Dilution Range	WB=1:500-2000, IHC-P=1:400-800, IF=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P analysis of mouse cerebellum using...



IHC-P analysis of rat cerebellum using M...



IHC-P analysis of human brain using MBP ...