

Product Datasheet

EZH2 antibody (orb783423)

Description	EZH2 antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	Recombinant human EZH2
Target	EZH2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1,IgK
Clonality	Monoclonal
MW	82 kDa
Uniprot ID	Q15910
Dilution Range	WB=1:500-2000, IHC-P=1:400-800, ICC=1:25, Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.