

Product Datasheet

LXR Alpha antibody (orb783426)

Description	LXR Alpha antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human LXR alpha (60-110/447 aa)
Target	LXR alpha
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	50 kDa
Uniprot ID	Q13133
Dilution Range	WB=1:500-1000, IHC-P=1:50-200, ICC=1:50, Flow-Cyt=1:50
Expiration Date	12 months from date of receipt.