



Product Datasheet

ANO1 antibody (orb783431)

Description	ANO1 antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant human ANO1/TMEM16A (37-333/986aa)
Target	ANO1/TMEM16A
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	114 kDa
Uniprot ID	Q5XXA6
Dilution Range	IHC-P=1:400-800
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human colon
carcinoma ...



IHC-P
analysis of
rat stomach
using ANO1...