



Product Datasheet

Lamin A antibody (orb783433)


Description	Lamin A antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human lamin A (1-100/664 aa)
Target	Lamin A/C
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	69/62 kDa
Uniprot ID	P02545
Dilution Range	WB=1:500-1000, IHC-P=1:50-500, ICC=1:50-100
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Pancreas (Rat); Lane ...



IHC-P analysis of mouse liver using Lami...



IHC-P analysis of human liver using Lami...