

Product Datasheet

Arg-vasopressin antibody (Biotin) (orb783435)

Description	Arg-vasopressin antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Arg-vasopressin(Cys-Tyr-Phe-Gln-Asn-Cys-Pro-Arg-Gly-NH ₂ , Disulfide bridge: 1-6)
Target	Arg-vasopressin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01185
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.