

Product Datasheet

Ubiquitin antibody (HRP) (orb783458)

Description	Ubiquitin antibody (HRP)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	ICC, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Ubiquitin
Target	Ubiquitin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P0CG48
Dilution Range	WB=1:500-1000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.