

Product Datasheet

CD33 antibody (FITC) (orb783466)

Description	Rabbit polyclonal antibody against CD33 conjugated to FITC
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from mouse CD33
Target	CD33
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40/25kDa
Uniprot ID	Q63994
Dilution Range	IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.