


Product Datasheet

Lamin A/C antibody (FITC) (orb7835)

Description	Mouse monoclonal antibody against Lamin A/C conjugated to FITC
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC
Immunogen	Recombinant bacterially expressed human lamin A
Target	Lamin A/C
Preservatives	Purified IgG conjugated to Fluorescein using Innova Biosciences Lightning-Link, supplied in Phosphate buffered saline (PBS) containing 0.09% Sodium azide
Concentration	1mg/ml
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P48678, P02545, P48679
Entrez	60374, 16905, 4000
Dilution Range	Optimal antibody dilution should be determined by titration
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of using Lam...