

Product Datasheet

A2793 (orb783559)

Description	A2793 is an efficient TWIK-related acid-sensitive K ⁺ channel (TASK)-1 inhibitor, with an IC ₅₀ of 6.8...
Note	For research use only
Formula	C ₁₃ H ₁₂ ClNO ₃
MW	265.69
CAS Number	[88349-90-0]
SMILES	O=C(OCC)COC1=C2N=CC=CC2=C(Cl)C=C1
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.