

Product Datasheet

CCR8 antibody (orb783576)

Description	Mouse monoclonal antibody to CCR8
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Recombinant Human CCR8 protein
Target	CCR8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	hIgG1
Clonality	Monoclonal
Uniprot ID	P51685
Expiration Date	12 months from date of receipt.