

Product Datasheet

Novel Coronavirus Spike glycoprotein(S) antibody (Biotin) (orb783581)

Description	Mouse monoclonal antibody against Novel Coronavirus Spike glycoprotein(S) conjugated to Biotin
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Novel Coronavirus Spike glycoprotein(S) (319-541aa)
Target	S
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	VHH fusion with human IgG1 Fc
Clonality	Monoclonal
Uniprot ID	P0DTC2
Dilution Range	Recommended dilution: ELISA:1:10000-1:50000
Expiration Date	12 months from date of receipt.