

Product Datasheet

Novel Coronavirus Spike glycoprotein(S) antibody (FITC) (orb783582)

Description	Mouse monoclonal antibody against Novel Coronavirus Spike glycoprotein(S) conjugated to FITC
Reactivity	Virus
Conjugation	FITC
Immunogen	Recombinant Human Novel Coronavirus Spike glycoprotein(S) (319-541aa)
Target	S
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	VHH fusion with human IgG1 Fc
Clonality	Monoclonal
Uniprot ID	P0DTC2
Expiration Date	12 months from date of receipt.