



Product Datasheet

PNKD antibody (orb78364)

Description	Rabbit polyclonal antibody to PNKD
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 164-342aa of human PNKD protein
Target	PNKD
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100µL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	43kDa
Dilution Range	WB: 1:500-1:2,000. IHC: 1:100-1:200. ELISA: 1:20,000-1:80,000
Expiration Date	12 months from date of receipt.



Western blot analysis of fetal liver lys...



Immunohistochemical analysis of formalin...