

Product Datasheet

ALPHA-ADR antibody (orb784460)

Description	ALPHA-ADR antibody
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Arabidopsis thaliana ALPHA-ADR protein
Target	ALPHA-ADR
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8LPL6
Expiration Date	12 months from date of receipt.