

Product Datasheet

3BETAHSD/D2 antibody (orb790050)

Description	3BETAHSD/D2 antibody
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Arabidopsis thaliana 3BETAHSD/D2 protein
Target	3BETAHSD/D2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q67ZE1
Expiration Date	12 months from date of receipt.