

Product Datasheet

MED19A antibody (orb792114)

Description	MED19A antibody
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Arabidopsis thaliana MED19A protein
Target	MED19A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9FMP0
Expiration Date	12 months from date of receipt.