

Product Datasheet

MAP70.3 antibody (orb792318)

Description	MAP70.3 antibody
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Arabidopsis thaliana MAP70.3 protein
Target	MAP70.3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9ZUA3
Expiration Date	12 months from date of receipt.