


Product Datasheet

CD71 antibody (orb79333)

Description	Goat polyclonal antibody to CD71
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Transferrin receptor 1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 90kDa band observed in Human Breast lysates (calculated MW of 84.9kDa according to NP_003225.2). Recommended concentration: 1-3 µg/ml. Experiment Notes: Reported variants represent identical protein: NP_003225.2, NP_001121620.1.
Protein Sequence	TRFSLARQVDGNSH
Clonality	Polyclonal
MW	84.9
Entrez	7037
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human breast ly...