


Product Datasheet

Myeloperoxidase antibody (orb79334)

Description	Goat polyclonal antibody to Myeloperoxidase
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	myeloperoxidase
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 85kDa band observed in Human Peripheral Blood Mononucleocytes lysates (calculated MW of 83.9kDa according to NP_000241.1). Recommended concentration: 0.3-1 µg/ml. Experiment Notes: The immunizing peptide represents the N terminus of the mature protein.
Protein Sequence	DKAYKERRESIQ
Clonality	Polyclonal
MW	83.9
Entrez	4353
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human Peripheral...