


Product Datasheet

CHRNE antibody (orb79344)

Description	Goat polyclonal antibody to CHRNE
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ACHRE / CHRNE
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 55kDa band observed in Mouse Fetal Skeletal Muscle lysates (calculated MW of 54.9kDa according to NP_033733.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	DQEATGEEVSD
Clonality	Polyclonal
MW	54.7
Entrez	1145
Dilution Range	ELISA: 1:1000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse Fetal skeletal muscle lysates.