


Product Datasheet

IL15 antibody (orb79353)

Description	Goat polyclonal to IL15
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	IL-15 (aa70-82)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoforms (NP_000576.1; NP_751915.1).WB: Approx 23kDa band observed in Human Bone Marrow lysates (calculated MW of 18.1kDa according to NP_000576.1).The observed molecular weight corresponds to the glycosylated form.Recommended concentration:1-3µg/ml.Experiment Notes: This antibody is expected to recognize both reported isoforms (NP_000576.1,NP_751915.1).
Protein Sequence	DATLYTESDVHPS
Clonality	Polyclonal
MW	18.1; 14.9
Entrez	3600
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human Bone Marrow lysates.