

Product Datasheet

PLA2-ALPHA antibody (orb793565)

Description	PLA2-ALPHA antibody
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Arabidopsis thaliana PLA2-ALPHA protein
Target	PLA2-ALPHA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8S8N6
Expiration Date	12 months from date of receipt.