


Product Datasheet

ADRA1B antibody (orb79357)

Description	Goat polyclonal antibody to ADRA1B
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ADRA1B (aa252-264)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 60kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 56.8kDa according to NP_000670.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	TLRIHSKNFHEDT
Clonality	Polyclonal
MW	56.8
Entrez	147
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human Cerebellu...