

Product Datasheet

NPF5.2 antibody (orb794589)

Description	NPF5.2 antibody
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Arabidopsis thaliana NPF5.2 protein
Target	NPF5.2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9FNL7
Expiration Date	12 months from date of receipt.