


Product Datasheet

Dopamine Receptor D1 antibody [Out of stock] (orb79896)

Description	Rabbit polyclonal antibody to Dopamine Receptor D1
Species/Host	Rabbit
Reactivity	Canine, Human, Monkey, Porcine
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic 15 amino acid peptide from 3rd cytoplasmic domain of human DRD1. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Monkey, Marmoset, Elephant, Panda, Dog, Bovine, Pig, Opossum (100%); Mouse, Hamster, Bat, Lizard (93%); Rat, Turkey, Zebra finch, Chicken (87%); Xenopus (80%).
Target	DRD1
Preservatives	PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Entrez	1812
Dilution Range	IHC-P: 8 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...