

Product Datasheet

Dopamine Receptor D1 antibody [Out of stock] (orb79897)

Description Rabbit polyclonal antibody to Dopamine Receptor D1

Species/Host Rabbit

Reactivity Human, Monkey

Conjugation Unconjugated

Tested Applications IHC-P

Immunogen Synthetic 17 amino acid peptide from C-Terminus of human DRD1. Percent identity with other species by BLAST analysis: Human, Gorilla, Monkey, Marmoset (100%); Gibbon (94%); Rat, Elephant, Pig (88%); Sheep, Hamster, Panda, Bovine (82%).

Target DRD1

Preservatives PBS, 0.1% Sodium Azide

Concentration 1 mg/ml

Storage Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.


Note For research use only

Clonality Polyclonal


Entrez **1812**

Dilution Range IHC-P: 5-10 µg/ml


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunohistochemical staining of Human Lu...



Immunohistochemical staining of human co...