

Product Datasheet

cyp-13A3 antibody (orb799413)

Description	cyp-13A3 antibody
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Caenorhabditis elegans cyp-13A3 protein
Target	cyp-13A3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q27517
Expiration Date	12 months from date of receipt.