

Product Datasheet

Y17G7B.18 antibody (orb799440)

Description	Y17G7B.18 antibody
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Caenorhabditis elegans Y17G7B.18 protein
Target	Y17G7B.18
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9U2R0
Expiration Date	12 months from date of receipt.