


Product Datasheet

Human beta Defensin 1 protein (orb80001)

Description	Beta Defensin-1 Human Recombinant produced in E.coli is a single, non-glycosylated polypeptide chain...
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	BD 1
Preservatives	The Human BD-1 was lyophilized from concentrated (1mg/ml) solution containing 20 mM PBS pH-7.4 and 130mM sodium chloride.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	GNFLTGLGHR SDHYNCVSSG GQCLYSACPI FTKIQGTCYR GKAKCKK.
Purity	> 98.0% as determined by RP-HPLC and analysis by SDS-PAGE
Solubility (25°C)	It is recommended to reconstitute the lyophilized Beta Defensin-1 in sterile H ₂ O not less than 100ug/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human beta Defensin...