


Product Datasheet

Rat FGF basic protein (orb80025)

Description	Rat FGF basic protein
Conjugation	Unconjugated
Tested Applications	FA, HPLC, SDS-PAGE
Target	Rat FGF basic
Preservatives	FGF-b was lyophilized from mg/ml solution after extensive dialysis against 20 mM phosphate buffer, pH 7.4 and 130mM NaCl.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	PALPEDGGGA FPPGHFKDPK RLYCKNGGFF LRIHPDGRVD GVREKSDPHV KLQLQAEERG VVSIKGVCAN RYLAMKEDGR LLASKCVTEE CFFFERLESN NYNTYRSRKY SSWYVALKRT GQYKLGSKTG PGQKAILFLP MSAKS.
Purity	> 97.0% as determined by RP-HPLC and analysis by SDS-PAGE
Solubility (25°C)	It is recommended to reconstitute the lyophilized Rat bFGF in sterile H ₂ O not less than 100ug/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Rat FGF
basic prote...