


Product Datasheet

Canine IL3 protein (orb80060)

Description	IL3 Canine Recombinant produced in E.coli is a single, non-glycosylated polypeptide chain containing...
Conjugation	Unconjugated
Tested Applications	FA, HPLC, SDS-PAGE
Target	k9 IL 3
Preservatives	Lyophilized from concentrated (1mg/ml) solution in water containing 1X PBS pH7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	RPFSTDLPKQ YFTMINEIME MLNKSPSPSE EPLDSNEKET LLEDTLRPN LDVFLNASSK FHKNGLLIWN NLKEFLPLLP TPTPRGEPIS IMENNWGDFQ RKLKKYLEAL DNFLNFKNKP.
Purity	> 97.0% as determined by: (a) Analysis by HPLC. (b) Analysis by SDS-PAGE.
Solubility (25°C)	It is recommended to reconstitute the lyophilized IL3 in sterile H ₂ O not less than 100ug/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Canine IL3 protein